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**Area of application:** Non precious alloy based on cobalt for removable dentures acc. DIN EN ISO 22674, type 5.

**Product description:** Omegalloy M is an alloy for denture bases is distinguished by an outstanding fluidity, which grants to fill even the thinnest details of the frame (0.5 mm in wax; 0.3 mm in metal). Moreover its molecular structure allows to obtain smooth, compact surfaces with little oxide formation. It is highly corrosion resistant and is free of toxic elements like beryllium, indium and gallium ( $\leq 0.1\%$ ).

#### Technical data (guidelines):

Proof stress 0,2 % (MPa):	650
Elongation at rupture (%):	5.0
Modulus of elasticity (GPa):	220
Tensile strength (MPa):	890
Vickers hardness VH 10:	350
Density (g/cm <sup>3</sup> ):	8.3
Melting range (°C):	1.280-1.360
Casting temperature (°C):	1.340

**Investing and casting:** Omegalloy M can be used with phosphate-bonded CoCr investments such as OmegaCast uni/rp. Prepare the moulds for casting according to the Omegatech DP ring less system. Preheating range is 850-1.000°C. Omegatech DP recommend 920°C. Heat soak at casting temperature for 30 minutes.

The alloy is cast with an induction machine when all the cylindrical ingots have melted together and a final shadow runs over the molten metal. Considerable experience is required when melting high-grade alloys automatically, and it is most important that the heating is set exactly.

Accurate heat setting is essential to protect the properties of the alloy. Light the flame and adjust the pressure in order to have a blue nucleus about 4 mm long. Place some metal ingots into the preheated crucible. Move the flame over the crucible rotating the tip, in order to evenly heat the alloy. As soon as the ingots are wholly melted, and after disappearing of their shadows, start centrifugal action. Casting Omegalloy M repeatedly is not recommended.

During re-use pay attention that the alloys are from the same batch. Mix old and new alloys in a ratio of 1:1.

**Preparing and polishing:** Cool down the muffle to room temperature, Sandblast the framework with aluminum oxide 120-150 µm and continue the working cycle according to your personal method. Remove carefully any small casting blows, go over all the corners and edges with a round bur and then smooth the surfaces with rubber abrasives and polish them.

**Soldering and welding:** Presolder using Cerabond Universal Solder (5g, REF 20009; 25g, REF 20010). For Laser weld-ing use Omegatech laser (Ø 0,25 mm, REF 24525; Ø 0,35 mm, REF 24535; Ø 0,50 mm, REF 24550) welding wire.

**Safety hints:** Metal dust is hazardous to health. For finishing and sandblasting use a suitable extraction system and / or face mask. The dental alloy can affect MRI results. Each supply is identified by a lot-number. Take note of this number on the patient's file in order to trace down the product.

**Contraindications, side effects and interactions with other dental alloys:** In cases of hypersensitivity (allergy) to the constituents of the alloy, discontinue its use. A Patch-Test is thus advisable. In individual cases, hyper-sensitivity reactions (allergies) and electrochemically induced local dysaesthesia have been reported, such as changes in taste and irritation of the oral mucosa. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys.

#### Product description

Co	Cr	Mo	Other constituents*
62,0 %	31,0 %	5,0 %	< 1,0%

\*Other constituents: Si, C, Mn

Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.

**Storage conditions:** No special storage conditions are required.

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SILADENT Dr. Böhme & Schöps GmbH kennzeichnet die Handelspackung mit: **CE** 0124
**Manufacturer:** SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE 38644 Goslar
SILADENT Dr. Böhme & Schöps GmbH packaging carries the mark: **CE** 0124