

# myerson FlexPress™

automated digital injection system

to proof  
6.17.2010



## Instruction Manual



*more to myerson than teeth*

[www.myersontooth.com](http://www.myersontooth.com)

|  |           |
|--|-----------|
| <b>Installation of the Heat Guard</b>  | <b>2</b>  |
| <b>Section 1 Components</b>  | <b>3</b>  |
| <b>Section 2 Storage &amp; Handling</b>  | <b>4</b>  |
| <b>Section 3 Basic Operation</b>   | <b>4</b>  |
| <b>Section 4 Editing User Profiles</b>   |           |
| <b>Section 5 Processing using the “Run Stored Profile” option (Recommended Option)</b> | <b>5</b>  |
| <b>Section 6 Processing using the “Rerun Last Profile”</b>                             |           |
| <b>Section 7 Myerson DuraFlex Partial Injection Procedure</b>                          | <b>7</b>  |
| <b>Section 8 Myerson DurAcetal Framework Injection Procedure</b>                       | <b>14</b> |
| <b>Section 9 Warranty &amp; Returns</b>  | <b>19</b> |
| <b>Section 10 Technical Data</b>   | <b>19</b> |
| <b>Section 11 Preventative Maintenance &amp; Cleaning</b>                              | <b>20</b> |
| <b>Section 12 Service &amp; Repairs</b>  | <b>20</b> |
| <b>Section 13 Declaration of Conformity</b>  | <b>19</b> |
| <b>Section 14 Important Safety Warnings</b>  | <b>21</b> |
| <b>Section 15 Notes</b>  | <b>22</b> |

## Installation of the HEAT GUARD

To prevent damage, the **myerson FlexPress™** has been shipped with the heat guard not installed. To install the heat guard, place the heat guard over the corner posts with the warning label UP. Install the (4) black acorn nuts finger tight. The acorn nuts are simply holding the heat guard down, so they do not need to be wrenched tight.

The heat guard can be removed periodically for cleaning. Just follow the above steps in reverse order to remove the heat guard.

Thank you for your purchase. Feel free to call myerson with any questions.

# section 1~components

The following is a list of components that are included in the **myerson Flexpress™** Digital Injection Systems:

| Description   | Professional Kit<br>120v - FP-0101<br>230v - FP-0103 | Advance Kit<br>120v - FP-0100<br>230v - FP-0102 |
|---|--|---|
| myerson <b>Flexpress™</b>   | X  | X   |
| (1) Injection Flask <b>BS-0104</b>  | X  | X   |
| (1) 6 Pack myerson <b>DuraFlex</b><br>Medium Cartridges                   | X  | X   |
| (1) 6 Pack myerson <b>DurAcetal</b><br>Medium Cartridges                  | X  | X   |
| (1) 6 Pack myerson <b>VisiClear</b><br>Medium Cartridges <b>VC-MD-6PK</b> | X  | X   |
| Power Cord  | X  | X   |
| Air hose  | X  | X   |
| (1) Thermoplastic Model Separator <b>MS-0200</b>                          | X  | X   |
| TMS Swabs (10) <b>MS-0202</b>   | X  | X   |
| TMS Trays (10) <b>MS-0201</b>   | X  | X   |
| Release Compound <b>BS-0105</b>   | X  | X   |
| Polishing Kit <b>FN-0200</b>  | X  |   |
| Bur Finishing Kit <b>AP-0224</b>  | X  |   |
| (1) Stone Removal Tool <b>SC-0100</b>                                     | X  |   |
| (1) 25lb. Universal Stone<br><b>UNISTONE-25LB</b>                         | X  |   |

## section 2 ~ storage & handling

The **myerson Flexpress™** Digital Injection System should be stored in a cool dry place. The unit should not be used in an environment exceeding 40°C (104°F). Keep all combustibles away from unit while in use. Never place anything on the top of the unit.

## section 3 ~ basic operation

Place the **myerson Flexpress™** Digital Injection System on a clean dry surface. Insert the power cord into the receptacle on the side of the unit and plug the other end into a grounded 120V or 230v outlet.

To turn on the unit, press the rocker switch located on the left side of the unit above the power cord receptacle. The control panel will illuminate briefly then turn back off. The **myerson Flexpress™** is now in Standby mode. To turn on the controller press the **ON/STBY** button (lower left button). Below is what you should see in the blue LCD display.



From the Main Menu you can select any of 5 options by using the **UP & DOWN** arrow keys. Below is the list of options & a description for each:

1. **Run Stored Profile-** **myerson,LLC** pre-programs the **myerson Flexpress™** with processing parameters for **myerson DurAcetal®**, **myerson DuraFlex** and **myerson VisiClear™**. By selecting the pre-programmed material, the Melt Temp, Melt Time and Hold Time are already set for you. This is the recommended option to use when processing cases using the **myerson Flexpress™**. There is also (3) "Spare" locations for the customer to save their own settings.
2. **Rerun Last Profile-** By selecting this option the **myerson Flexpress™** will run the last set profile. This option is great when processing several cases of the same material.
3. **Edit User Profiles-** This option is where you can create your own processing profiles and save them for future use. The **myerson FlexPress** is capable of storing up to 7 unique processing profiles.
4. **Edit User Settings-** This option allows you to change settings such as temperature units (C/°F), turn off/on alarms, languages, etc.
6. **Configure from PC-** This option is used by repair personnel in the event of a warranty repair of the **myerson FlexPress**.

## section 4 ~ editing user profiles

The **myerson FlexPress** allows the user to save up to 7 custom processing profiles for easy switching between materials. The following procedure outlines how to create and save profiles in the **myerson FlexPress™** Digital Injection System. Note: **myerson, LLC** does not provide process settings for materials other than **myerson, LLC** supplied materials. Refer to the manufacturer supplied information for processing settings.

1. From the main menu, use the **UP/DOWN** buttons until **“Edit User Profile”** is displayed. Press the **ENTER** button to select.
2. Using the **UP/DOWN** buttons until the “Spare” location to be edited is displayed. Press the **ENTER** button to select.
3. The first parameter to be modified is the “Name” of the spare profile. Press the **ENTER** button and a flashing cursor should appear. Using the **LEFT/RIGHT** buttons to move the cursor to the character location you wish to modify. Once the cursor is in position, use the **UP/DOWN** buttons to scroll through the available characters. Press the **ENTER** button when finished.
4. Press the **DOWN** button and the “Melt Temperature” parameter will be displayed. Press the **ENTER** button and a flashing cursor will appear. Use the **UP/DOWN** buttons to change the Melt temperature to the desired value. Once finished press the **ENTER** button.
5. Repeat step 4 for the **“Melt Time”** and “Hold Time” parameters.
6. Once satisfied with the changes, press the **LEFT** button. The screen should display **“Saving Changes”**. The profile you just created will now appear in the **“Run Stored Profile”** option in the main menu.

## section 5 ~ Processing using the “Run Stored Profile”

### (Recommended Option)

1. From the main menu use the **UP or Down** buttons until the “Run Stored Profile” is displayed then press the enter button.
2. Next use the **UP/DOWN** buttons until profile you wish to run is displayed.
3. After selecting the profile the display will show the profile name. Pressing **ENTER** will scroll between profile name, melt temperature, melt time and hold time. It is recommended to check the setting before starting a profile.
4. Press the **START** button and the machine will begin running the profile. If the machine is at room temperature, it will take approximately 15 minutes to reach the processing temperature.
5. Prepare the cartridge by scoring a “+” using a utility knife. See the illustration to the right. This will create an even “Burst” of the cartridge.
6. Next, coat the sides of the cartridge with a light coat of release compound. This will prevent the cartridge from sticking to the walls of the heating chamber.
7. When the **myerson Flexpress** has calibrated itself at the Melt Temp an alarm will sound and the upper display will read “**LOAD**” signaling you to load the cartridge. At this time insert the cartridge into the heating chamber, center the flask and tighten using the flask strap and knobs. If you are not ready to start the profile at this time and want to turn off the alarm, press the UP button. This will silence the alarm, but the **myerson Flexpress™** will remain calibrated until you press **START** or put the unit in standby by pressing the STBY button.
8. Press the **START** button to begin the profile. Now let the **myerson Flexpress™** do the rest. The unit will begin counting down the Melt Time and once it expires the unit will automatically inject the material into the flask. The piston will stay in the up position for the duration of the Hold time then release.
9. The upper display will read “**DONE**”.
10. You can now remove the flask strap and knobs.
11. Grasp the flask using heat resistant gloves and press the **EXTEND/RETRACT** button. This will extend the piston and push out the remaining cartridge attached to the flask.



12. Press the **EXTEND/RETRACT** button again to retract the piston.
13. The **myerson Flexpress™** is now ready to process another case. Repeat steps 1-9 for your next injection.

## section 6 ~ processing using the “Re-run Last Profile”

1. From the main menu use the **UP or DOWN** buttons until the cursor is to the left of “Rerun Last Profile” then press the enter button.
2. The display will show the profile name and profile settings of the last injection. If this is not the profile you wish to run, press the LEFT arrow button to return to the main menu.
3. The **myerson Flexpress™** is now ready to process your case. See steps 5-13 of Section 5 for a complete injection.

## section 7 ~ myerson DuraFlex partial injection procedure

The following procedure outlines the steps necessary to create **myerson DuraFlex** Partials with the **myerson Flexpress™** Digital Injection System.

1. Evaluate the model for major undercuts and block out the undercuts. Due to the flexibility and adjustability of the material, very little if any blocking out will be needed.



2. Paint the model with myerson Thermoplastic Model Separator **MS-0200**.



3. Wax and set up the teeth just as you would for an acrylic partial.



4. Add the clasps to the partial. **myerson,LLC** recommends using 12 gauge sprue wax for creating clasps.







5. Coat the inside of the flask with a petroleum jelly. This will prevent the stone from adhering to the flask.
6. Invest the partial using regular Buff stone. When stone sets, place sprues using 3/16" (4.75mm) utility rope wax. Extend the sprue to the injection port in the flask.



7. Coat the stone with stone separator and bolt the top half of the flask in place. The utility wax should fill the injection port of the flask. This will keep the stone from leaking out during the second half of the investing process.



8. Pour up the flask using **myerson, LLC's** Investing Stone **INVSTONE-25LB**. This stone has been designed to prevent the teeth from shifting during injection. Notice the stone is not protruding from the widow in the flask.



9. When the last half of the stone is set, remove the bolts but leave the flask together and place in boiling water.
10. Once the wax has been removed, remove and place holes in the denture teeth. The holes in the teeth are very important. They create the mechanical retention necessary to lock the teeth in place as well as allow air to escape as the **myerson DuraFlex** is injected. Once the holes have been created, the teeth need to be placed back into the proper locations. A small dot of Super Glue on the incisal area of the teeth will hold them in place as the material is injected. Bolt the two halves together.




11. Next, turn on the **myerson Flexpress™** by pressing the **ON/STBY** button. Select the method of processing referred to in sections 4-6.
  
12. When the **myerson Flexpress™** has reached the preset temperature, a buzzer will sound. At this time, insert the desired **myerson DuraFlex** cartridge into the heating chamber with the crimped end down. A light coat of release compound BS-0105 on the cartridge will help eliminate the cartridge getting stuck in the heating chamber. Press the start button this will start the injection cycle.



13. Place the flask onto the **myerson Flexpress™**, centering the sprue hole of the flask over the center of the heating chamber. Lock the flask in place using the flask strap and knobs.



14. When the timer expires, the **myerson Flexpress™** will automatically inject the **myerson DuraFlex** cartridge. The piston will stay in the up position according to the hold time. Once the piston retracts, it is safe to unscrew the hold down plate and remove the flask.
15. To remove the flask, loosen the hold down knobs and remove the flask strap. Using heat resistant gloves, grasp the flask and press the **EXTEND/RETRACT** button (lower right button). This will extend the piston and push the crushed cartridge and flask out.
16. Press the **EXTEND/RETRACT** button again to retract the piston. The **myerson Flexpress™** is now ready for another injection.
17. After the appliance is divested, it can be placed in a stone remover to eliminate any extra stone that is adhering to the appliance. Notice how the Thermoplastic Model Separator **MS-0200** has created a smooth, clean surface on the underside of the appliance.



18. This is the finished **myerson DuraFlex** case. Ready for try in.



# Section 8 ~ myerson DurAcetal framework injection procedure

The following procedure outlines the steps necessary to create **myerson DurAcetal®** frameworks with the **myerson Flexpress™** Digital Injection System.

1. Block out the saddles on the model. This is in preparation for the acrylic that will be placed later. The thickness used is 1 to 1.5 mm thick. You will notice the large tissue stop cut out of the wax. Stability is the reason for the large tissue stop.



2. Below is the duplicated model. This model has been duplicated using **myerson** Universal Stone **UNISTONE-25LB**. The expansion stone accounts for the shrinkage that occurs when **myerson DurAcetal** cools.





3. Place myerson's Thermoplastic Model Separator **MS-0200** on the duplicated model.



**myerson**  
**Cosmetic Dental Materials**



4. Wax up the framework.



5. Coat the inside of the flask with a petroleum jelly. This will prevent the stone from adhering to the flask.
6. Invest the framework using regular Buff stone. When stone sets, place sprues using 3/16" (4.75mm) utility wax. Extend the sprue to the injection port in the flask.



7. Coat the stone with stone separator and bolt the top half of the flask in place. The utility wax should fill the injection port of the flask. This will keep the stone from leaking out during the second half of the investing process.
8. Pour up the flask using **myerson, LLC's** Investing Stone **INVSTONE-25LB**. Notice the stone is not protruding from the widow in the flask.




9. When the last half of the stone is set, remove the bolts but leave the flask together and place in boiling water.
10. Once all the wax has been removed, bolt the two halves together.





11. Next, turn on the **myerson Flexpress™** by pressing the **ON/STBY** button. Select the method of processing referred to in sections 4-6.
12. After the appliance is divested, it can be placed in a stone remover to eliminate any extra stone that is adhering to the appliance. Notice how the Thermoplastic Model Separator **MS-0200** has created a smooth, clean surface on the underside of the appliance.



13. Here you have your finished **myerson DurAcetal** appliance. The holes in the saddle areas were created using a standard round bur. This will allow the acrylic to mechanically retain to the **myerson DurAcetal** framework.



## Section 9 ~ myerson warranty & returns

**myerson** offers a three (3) year limited warranty for US and Canada against defects in material and workmanship. If you would like to submit a product for warranty return, please contact a **myerson** distributor or **myerson**. Include only the defective product submitted for warranty replacement; also include a letter describing in detail the problem, date of purchase of the product, as well as a return address and contact name. **myerson, LLC** is not responsible for return shipment costs. **myerson** will refurbish or replace product found to be defective of material or workmanship. If after the three (3) year warranty you are experiencing a problem with a product, please contact a **myerson** distributor or myerson.

## Section 10~ technical data

- Power Supply            120VAC 50/60Hz or 230VAC 50/60Hz
- Max Melting Temp    608F (320C)
- Dimensions            H: 17in, W: 11in, D: 10in (43cm,28cm, 25cm)
- Weight                 29lbs (13.15 kg)
- Air Pressure            100 psi Recommended; 150 psi Maximum
- Fuse                     10 amp, (5x20mm)

## Section 11 ~ preventive maintenance & cleaning

Use only approved materials in the **myerson Flexpress™** Digital Injection System. The use of other materials may cause the unit to fail and will void the manufactures warranty. See Section 10 for warranty information.

Keep heating chamber and flask area free from dirt and other foreign debris to prolong the life of the unit. Periodically clean the heating chamber with the cleaning brush. The rest of the unit can be cleaned with a damp rag and mild soap solution. Turn off and unplug the unit before cleaning. Always keep the unit in an upright position.

## Section 12 ~ Service & Repairs

For service and repair of the myerson **Flexpress™** Digital Injection

System, please call **888.401.1126** or go to [cdmdental.com](http://cdmdental.com).

For technical support contact **Kris Schermerhorn** CDT at **800.874.1047** or email at **Kfuldek@msn.com**

## Section 13 ~ Declaration of Conformity

Use only approved materials in the **myerson Flexpress™** Digital Injection System. The use of other materials may cause the unit to fail and will void the manufactures warranty. See Section 10 for warranty information.

Keep heating chamber and flask area free from dirt and other foreign debris to prolong the life of the unit. Periodically clean the heating chamber with the cleaning brush. The rest of the unit can be cleaned with a damp rag and mild soap solution. Turn off and unplug the unit before cleaning.

**Always keep the unit in an upright position.**

## Section 14 ~ important safety warnings

When using the **myerson Flexpress™** Digital Injection System, always follow basic safety precautions to reduce the risk of accidental injury due to burns or electrical shock.

Wear long sleeves, long pants, heat-resistant gloves, and safety goggles when using the unit. Do not expose heat-resistant gloves to hot metal areas for extended periods of time.

Keep out of the reach of children. Do not place unit where cord can be tripped over. Keep unit away from flammable material. It is normal for the unit to become quite warm during use.

**Do not leave unit unattended with the power on.**

The **myerson Flexpress™** Digital Injection System is for indoor use only. Do not use outdoors or in wet or damp locations. Never immerse unit in water or other liquids. Use only on a clean, dry work surface. Excessive dust and dirt will reduce the life of the unit.

Use only recommended accessories and materials with this unit. Using materials other than those recommended by **myerson, LLC** may result in accidental injury and permanent damage to the unit. When not in use, unit should be stored in a dry secure location out of the reach of children.

Manufacturer assumes no responsibility for consequential or indirect damages from the use of this product.

Use this product only as specified in this manual. Using this device in a manner not specified may result in injury or hazards.

**WARNING** Do not turn off power when cooling fan is running.

# section 15 ~ notes

**myerson**  
the tooth company

more to myerson than teeth  
www.myersontooth.com

**myerson SPECIAL**

Premium quality teeth hand-crafted to be virtually indistinguishable from natural teeth



**myerson DB PLUS**

Time proven formula reinforced with sub-micron silica to enhance hardness




**KENSON**

"The Professional Choice" for over 50 years; Life-like aesthetics, unparalleled value




**myerson-ema**  
...Sleep Better




ema® is a trademark of frantz design inc.

**myerson.**  
Cosmetic Dental Materials

**myerson DuraFlex**  
flexible partial dentures



**myerson DurAcetal**  
durable cosmetic appliances



**myerson VisiClear**  
clear clasps & frameworks



You could be using the **myerson ema oral appliance** for the treatment of sleep related breathing disorders. **ema's** patented design uses interchangeable elastic straps attached to pressure-formed custom upper and lower trays, which allow the gradual advancement of the mandible until treatment is successful. **myerson ema oral appliance** is available to medical and dental professionals through authorized dental laboratories.

**jim please check this copy....**



*more to myerson than teeth*  
[www.myersontooth.com](http://www.myersontooth.com)

**MANUFACTURED BY**

myerson company limited  
3 trinity avenue  
laventille, trinidad & tobago  
telephone 868 623 1007  
facsimile 868 627 4594



myerson LLC  
311 north aberdeen  
chicago, IL 60607-1203 USA  
toll free phone 800 423 2683  
toll free fax 800 424 2928  
telephone 312 432 8200  
facsimile 312 563 9535

**EUROPEAN REPRESENTATIVE**

myerson limited  
8 crystal way, harrow middx,  
HA1 2HP, united kingdom  
telephone +44 (0) 20 8863 9044  
facsimile +44 (0) 20 8861 3091